

## Selection® Ulnar deviation Splint

Improve finger alignment to enhance  
hand function and functional grip



## Selection® Ulnar deviation splint

Selection® Ulnar Deviation Splint is intended to improve finger alignment and enhance hand function and functional grip. The fabric can be trimmed without fraying, and therefore allows further customization if required.

The easy fastening surface allows the hook to be placed anywhere on the splint to optimise the directional pull needed. The orthosis is designed so that the material extends far to the ulnar side at 5. MCP to support the best possible direction of pull and distribute the pressure.

The Ulnar Deviation splint supports the arches of the hand whilst providing re-alignment, stabilisation and can contribute to the management of pain. Should you need to provide wrist stabilisation for your patient, this splint can also be worn in combination with a Selection wrist splint.

### Intended use

For patients with inflammatory and degenerative diseases, such as rheumatoid arthritis, that has resulted in ulnar deviation. Selection Ulnar Deviation can also be used for hypermobile joints and weak muscles to provide stability. Not intended for use with fixed contractures in the hand.

Art. No. 35218



Scan or click for sizes and measurements!

## Improve finger alignment to enhance hand function and functional grip

This hand splint is exclusively designed and produced using our unique three layer Selection fabric. With its low profile design, it meets the contours of the hand without impeding function.

- Meets the shape and contours of the hand
- Breathable fabric
- Fabric can be trimmed
- Easy fastening
- Can be worn with Selection Wrist Splint
- Machine washable

## Developed for optimal function and comfort

### Selection® - a full range of hand and wrist orthoses

Our Selection range is created for anatomically compliant orthoses that provide support and stability where needed, but at the same time are flexible and allow movement so that the body part can be functional. Our mission with Selection is for the products to feel like a second skin and to be customisable for a perfect fit for everyone.


The soft interior and edges, together with the breathable and quick drying material, provides a comfortable and cool fit. Selection orthoses has lacing for customized size and support and can be cut for further customization without fraying. The low profile of the products allows them to be used under clothes and gloves.

Tested to be used closed to the body according to SGS.



### The Selection Material

All our materials are carefully tested at internationally approved testing institutes to ensure that they contain no substances that can be harmful to either adults or children. The Velcro straps used in the products are carefully selected to suit their particular function and has also been tested to ensure long lasting durability.



### Selection Powernet

The function carrier, the inner textile mesh structure, gives the material its perfect elasticity and recovery for good consistent support.

### Selection Air

Thin, soft and 100% air-permeable material. Used in selected parts of the orthosis where the textile should feel invisible and not retain heat. Gives extra smooth function.

### An excellent combination

The combination provides an anatomical orthosis that allows unrestricted movement and stability. A perfect fit that is soft inside and durable outside. Elegant with high comfort.



## **Support for Better Life**

Everyone should be able to live their life to the fullest, regardless of their mobility challenges. With innovative solutions developed in close collaboration with healthcare professionals and patients, we strive to provide Support for Better Life.

# allard

Allard International

Karbingatan 38, 254 67 Helsingborg, Sweden | +46 42 25 27 00 | [info@allardint.com](mailto:info@allardint.com) | [allardint.com](http://allardint.com)