



SELECTION[®] RADIALIS FINGER & THUMB PART


The SELECTION[®] range presents two new products, Radialis Finger and Radialis Thumb. These are sold as separate items to make it possible to combine them optimally. They are intended to be used when a Radialis Paresis is present (drop hand caused by damage to the radialis nerve) to assist functional gripping ability whilst the nerve recovers. These are intended to be combined with SELECTION[®] Wrist orthosis model numbers 35204, 35205 and 35209.



The Thumb Part, shown above, is attached to SELECTION[®] Wrist Soft 35204.

Item no.	Size	Side	Measure
35216 1011	Small	Left	3 - 4 cm
35216 1013	Large	Left	5 - 6 cm
35216 2011	Small	Right	3 - 4 cm
35216 2013	Large	Right	5 - 6 cm

Measure: Length joint centre MCP1 – joint centre IP1.





SELECTION®

RADIALIS

FINGER & THUMB PART

Radialis Finger Part comprises two sections that makes it easier to put on and adapt to the best fit. The Finger Part is attached to the soft plush material on the dorsal side of the orthosis. Simple to attach and comfortable to use. Can preferably be combined with Radialis Thumb Part.



The Finger Part, shown above, is attached on a SELECTION® Wrist Soft 35204. The image to the left shows the Thumb and Finger Parts attached to the orthosis.

Item no.	Size	Side	Measure
35217 1011	Small	Left	17 - 20 cm
35217 1012	Medium	Left	20 - 23 cm
35217 1013	Large	Left	23 - 27 cm
35217 2011	Small	Right	17 - 20 cm
35217 2012	Medium	Right	20 - 23 cm
35217 2013	Large	Right	23 - 27 cm

Measure: Circumference MCP 2-5

