

SELECTION®

Soft Knee Bracing in a totally new way



INTENDED USE

These products are intended for patients suffering from mild knee pain, who may benefit from a light degree of compression and some warmth. These products are only intended for one person use.

INDICATIONS

Knee pain in weight bearing, minor injuries such as stress in ligaments or mild sprains.

CONTRAINDICATIONS

Severe knee injuries such as ACL/PCL ruptures.
Caution: If severe swelling or blood circulation problems exist.

INTENDED USE

This product is intended for patella support of various kinds, as the patella pad in the brace can be moved into different positions. When positioning the patella pad on the side, the brace can be used by patients suffering from lux/subluxation of patella. When positioning the patella pad below the patella, it can be used by patients suffering from patella tendon pain, Osgood Schlatters, Jumpers knee or Patellafemoral Pain Syndrome.

INDICATIONS

Patella luxation (hyper mobility), patella tendon pain, jumpers knee, Osgood Schlatter, Patellafemoral Pain Syndrome.

CONTRAINDICATIONS

None known. Caution: If severe swelling or blood circulation problems exist.



INTENDED USE

These products are intended for patients suffering from various knee pain, as the products contain many different support structures. They give assistance in reducing hyperextension, adding M-L stability, and can assist in A-P instability. Utilizes the ROM-hinge for post surgical treatments.

INDICATIONS

Mild to moderate knee hyperextension, weak ACL, partial ACL rupture, mild to moderate M-L instability, mild osteoarthritis, post surgery.

CONTRAINDICATIONS

Contact sports or other activities where the knee is exposed to a high degree of damage and the risk for ACL/PCL injuries is severe.

Caution: If severe swelling or blood circulation problems exist.

allard^{INT}

SUPPORT FOR BETTER LIFE



Minor



Optima




Patella

SELECTION®

To be able to give as good fit as possible the sizing chart for the SELECTION range of products has 3 different measuring points - see table.

The two WRAP-AROUND-MODELS also have an additional PLUS size to accommodate larger thigh circumference.

SELECTION® Minor
SELECTION® Optima
SELECTION® Patella



Minor	Optima	Patella	Size	15 cm below Patella	Mid Patella	15 cm above Patella
35700 0010	35702 0010	35704 0010	X-Small	28-31 cm	30-33 cm	34-39 cm
35700 0011	35702 0011	35704 0011	Small	31-34 cm	33-36 cm	39-43 cm
35700 0012	35702 0012	35704 0012	Medium	34-37 cm	36-39 cm	43-47 cm
35700 0013	35702 0013	35704 0013	Large	37-40 cm	39-42 cm	47-51 cm
35700 0014	35702 0014	35704 0014	X-Large	40-43 cm	42-46 cm	51-56 cm




Minor Wrap



Optima Wrap

SELECTION® Minor Wrap
SELECTION® Optima Wrap



Minor Wrap	Optima Wrap	Size	15 cm below Patella	Mid Patella	15 cm above Patella
35701 0010	35703 0010	X-Small	25-28 cm	27-30 cm	31-36 cm
35701 0011	35703 0011	Small	28-31 cm	30-33 cm	36-40 cm
35701 0115	35703 0115	Small plus	28-31 cm	30-33 cm	40-44 cm
35701 0012	35703 0012	Medium	31-34 cm	33-36 cm	40-44 cm
35701 0125	35703 0125	Medium plus	31-34 cm	33-36 cm	4-48 cm
35701 0013	35703 0013	Large	34-37 cm	36-39 cm	44-48 cm
35701 0135	35703 0135	Large plus	34-37 cm	36-39 cm	48-53 cm
35701 0014	35703 0014	X-Large	37-40 cm	39-43 cm	48-53 cm
35701 0145	35703 0145	X-Large plus	37-40 cm	39-43 cm	53-60 cm