



SELECTION[®] **MULTIOPEN**

DOUBLE locking function

SELECTION[®] Multiopen is a semi-stable wrist support with double locking function. It can be opened up completely and secured loosely with Velcro and then pulled closed firmly with the lacing system.

The patient can choose to either open it completely each time the orthosis is taken on/off or only use the lacing function.

The orthosis is low cut to follow the hand's natural bending folds and the recess around the thumb allows the the muscles to work freely.

The two-part lacing and locking in combination with Velcro enables adjustment of the width for a wider forearm.

As the lacing is positioned on the ulnar side, the risk of the thumb being positioned incorrectly decreases. The orthosis has an adjustable thumb strap that attaches from the volar side and fastens on the dorsal side.

Through the construction of the material, it is possible to cut, using scissors, to shape for optimal fit. Made in an airy and comfortable SELECTION[®] material with all the functions that this entails.

The orthosis is 22cm in length, i.e. 4cm longer than a SELECTION[®] Soft 35204.

It has a malleable aluminum stay on the volar side which is removable.

SELECTION®

MULTIOPEN

When Multiopen is applied the first time it should be opened out completely. You will see the opening in the top of the adjacent picture. Attach toward the velour loosely.



The Velcro does not need to match the velour surface exactly.

Here the lacing has been tightened. These two locking functions together gives an orthosis that also adapts according to the width of the arm.



Item No. Left	Item No. Right	Size	Measure
35209 1010	35209 2010	X-Small	< 11 cm
35209 1011	35209 2011	Small	11 - 14 cm
35209 1012	35209 2012	Medium	14 - 17 cm
35209 1013	35209 2013	Large	17 - 20 cm
35209 1014	35209 2014	X-Large	20 - 23 cm
35209 1015	35209 2015	XX-Large	23 - 26 cm

